



# Nano

## Portable Oxygen Concentrator



### Model Part Number

**855**

100 to 240 VAC  
50 to 60Hz

---

### Specifications

**Average Oxygen Content** 90% (+6% / -3%)  
at all settings

**Rechargeable Battery** 14.4 VDC Voltage /  
6.8Ah Rated Capacity

**Breathing Frequency** 10 to 40 BPM

**Weight** 4.7 lbs / 2.1 kg

**Dimensions** 8.3"W x 3.5"D x 6.5"H  
23 cm x 9 cm x 17 cm

2010-8100 REV A

**Light Weight!**

**Take it Anywhere!**

**Bright LCD Screen for Easy Viewing!**

**Easy Touch Button Controls!**

**Convenient Carrying Bag Included!**

**NIDEK**  
*Medical*



MADE IN  
**USA**

**CE**  
0413